

Assembly of Image Guided Laptop Stands for KneeTrac CPTs,

Purpose

This document is a production procedure for the Knee Trac System to make sure that the system will be assembled as released.

The assembly has to be carried out by STRYKER LEIBINGER engineers only !

Identification

Project Name / - no.	KneeTrac 1.0 / #1647
End User / Hospital	
Option (Duracon or Scorpio)	
Part Number	
Serial Number	
Responsible for assembly (Print Name)	

Appendices

Distributed to

☒ DHR☐ Other _____

Summary

- ☐ OK The system is complete and was strictly assembled according to the Assembly Instructions. (page 4) All OK-fields were checked.
- ☐ OK (minor deviations) (deviations were listed on last page). These deviations wont cause any functional limitations or any hazards !
- ☐ NOT OK correct and complete assembly is not possible (Reasons were listed on last page).

Approvals
(Responsible
Engineer)

(Print Name)

(Signature)

(Date)

Continued on next page

100701-4254/650

Assembly of Image Guided Laptop Stands for KneeTrac CPTs, Continued

Part List

Collect serial numbers and check completeness and appearance					
No	Description	Qty	Serial no / ID	OK ?	
				No	Yes
1	Laptop Stand w/ Modified Laptop tray	1		<input type="checkbox"/>	<input type="checkbox"/>
2	6000-500-28 IGT Localizer (includes power cord and null-modem cable)	1		<input type="checkbox"/>	<input type="checkbox"/>
3	6000-500-006 IGT Sensor Array (includes camera cable)	1		<input type="checkbox"/>	<input type="checkbox"/>
4	Dell Laptop CPU	1		<input type="checkbox"/>	<input type="checkbox"/>
5	Dell Port Replicator (includes power cord)	1		<input type="checkbox"/>	<input type="checkbox"/>
6	LCD Monitor (includes power cord, VGA cable)	1		<input type="checkbox"/>	<input type="checkbox"/>
7	Wiremold Power Strip (includes mounting brackets)	1		<input type="checkbox"/>	<input type="checkbox"/>
8	Isolation Transformer (includes mounting brackets and 230 V fuses)	1		<input type="checkbox"/>	<input type="checkbox"/>
9	Isolation Transformer Mounting Plate	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
10	6000-115-5 Universal Mount	2	N/A	<input type="checkbox"/>	<input type="checkbox"/>
11	6000-500-27 Mounting Plate (Sensor array)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
12	6000-102-27 Mounting Plate (Monitor)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
13	6000-115-155 Arm Assembly (Includes cable housing cover, 8 screws)	2	N/A	<input type="checkbox"/>	<input type="checkbox"/>
14	6000-115-157 Arm Channel (Monitor)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
15	6000-115-216 Arm Channel (Sensor array)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
16	Arm Bushing	2	N/A	<input type="checkbox"/>	<input type="checkbox"/>
17	0001-143-000 Screws (to fasten universal mount to mounting plate)	8	N/A	<input type="checkbox"/>	<input type="checkbox"/>

Continued on next page

100101-4254/660

Assembly of Image Guided Laptop Stands for KneeTrac CPTs, Continued

Part List (continued)

Nr	Description	#	Serial no / ID	OK ?	
				No	Yes
18	0002-106-000 Screws (to fasten monitor mounting assembly to monitor)	4	N/A	<input type="checkbox"/>	<input type="checkbox"/>
19	Reverse IEC320 adapter kit (for power strip to isolation transformer)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
20	Image Guided Cart power cord	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
21	European plug adapter kit (for power cord to wall)	1	N/A	<input type="checkbox"/>	<input type="checkbox"/>
22	Ground wire w/ Ring Terminal Connectors	1		<input type="checkbox"/>	<input type="checkbox"/>
23	#10-32 x 1/2" Machine Screws	12		<input type="checkbox"/>	<input type="checkbox"/>
24	#10-32 Nuts	4		<input type="checkbox"/>	<input type="checkbox"/>
25	#10 Star Washers	21		<input type="checkbox"/>	<input type="checkbox"/>
Software					
26	GhostVersion:	1		<input type="checkbox"/>	<input type="checkbox"/>
27	KneeTrac -SW:	1		<input type="checkbox"/>	<input type="checkbox"/>
28	Printed User's Manual V1.0	1		<input type="checkbox"/>	<input type="checkbox"/>

Smart Tools & Navigational Hardware					
No	Description		Qty	OK ?	
				No	Yes
1.	6003-011-000 KneeTrac Pointer	ID-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>
2.	6003-005-000 Tibial / Pelvic Tracker	ID-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>
3.	6003-010-000 Femoral Tracker	ID-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>
4.	6003-025-000 Resection Plane Probe	PT-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>
5.	6003-019-000 Trocar for Pins	PT-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>
6.	6003-020-000 Sleeve with Handpiece	PT-no.: _____	1	<input type="checkbox"/>	<input type="checkbox"/>

FOOT-4254/660